

Build a sentence

Step-by-step instructions for Build a sentence in Express Me.

Watch the video, follow the steps, or open the complete transcript for a text version of the instructions.



Scan for video

Step-by-step instructions

1

Enable sentence building

Open **Settings** and enable the option to build a sentence with selected words.

2

Add tiles to the sentence

Select a tile to add it to the sentence. The sentence is built in the section at the bottom of the tile page.

3

Clear the sentence

Click the **Delete** button to clear the sentence.

Complete transcript

▼ Read the complete video transcript

I'll show you how to build a sentence with words in the app. To do this, first open settings, and make sure the option to build a sentence with selected words is enabled.

When you've done this, you'll see a new section at the bottom of the tile page where the sentence is being constructed. Now when you select a tile, it will be added to the sentence.

By pressing the button to delete the sentence, the sentence will be cleared.